

STG – SIP to TAP Gateway

• Rauland® SIP Protocol Monitoring • TAP Protocol Output •



“Enhancing Mobile Alert Response”

Visit Our Website at:

www.WirelessMessaging.com

2630 National Dr. Garland TX. 75041 Ph: 1.800.373.1466

The WaveWare **STG** SIP to TAP Gateway allows you to Provide TAP Protocol Outputs from your Rauland™ Nurse Call System using a Single Embedded System

STG Supports up to One Ethernet Based TAP Output and up to Four RS-232 TAP Outputs (Requires USB to Serial Adapters)

Features:

The Superior Features of the WaveWare **STG** SIP to TAP Gateway Include:

- **Ethernet** Port – Supports Rauland™ SIP Monitoring and TAP Output
- **Serial** Ports – Qty. up to 4 TAP outputs using USB to Serial Adapters
- Simultaneous TAP Output via One or More USB Ports and One Ethernet Socket Connection
- Logging of input and output transactions
- Supports TAP 1.8 and Rauland SIP Interface Protocols
- Configurable using Web Browser
- System Status and Logging via Web Browser
- Discoverable on LAN using WaveWare Discover and Reset App
- Programmable IP based and Serial Com Port Settings and Protocols



P/N= **0250-1200**

System Specifications:

- Model Numbers: **STG SIP to TAP Gateway**
 - Part Number: **0250-1200**
 - **Ethernet** Port: RJ-45 Female
 - 4 USB Ports (requires Belkin USB-to-Serial Portable Adapters)
 - Interface Protocols: Rauland SIP, TAP over Ethernet, TAP over Serial Port
 - Configuration using Web Browser
 - Power Supply: 19V 65W Wall Mount AC/DC Universal Switching Adaptor
 - Dimensions: Approx: 2” (H) x 5” (W) x5” (D)
- Weight: Approx: 2.5 lbs

Packaging:

The **STG** SIP to TAP Gateway Comes Packaged with Everything you need to Configure and Install the System. Including:

- Qty. One – STG Server on an Intel NUC computer running Linux, with AC/DC Power Adapter
- Qty. One – Belkin USB to RS-232 Serial Adapter
- Qty. One - Standard Ethernet Cable